

TISSUE RUN UP FOR LR WHITE RESIN

Gloves must be worn throughout the procedure. All work should be done in a fume hood.

- Fix tissue in 4% paraformaldehyde, and 0.1% glutaraldehyde in 1x PBS (pH 7.4) for 4-6 h..
- Wash in 3 changes 1x PBS for 15 minutes each.
- Dehydrate in graded series of alcohol for 15 minutes each:
 - 30% ethanol
 - 50% ethanol
 - 70% ethanol
 - 90% ethanol
 - 100% ethanol x 3
- Infiltrate with LR White resin in 2 one hour changes.
- Embed in LR White resin.

IMMUNOGOLD LABELING OF LR WHITE EMBEDDED SECTIONS

- Block with normal goat serum (1:20) for 30 minutes in a humid chamber at room temperature.
- Incubate with primary antibody for 60 minutes in a humid chamber at room temperature.
- Rinse in BSA three times for 5 minutes each.
- Incubate with secondary antibody for 60 minutes in a humid chamber at room temperature.
- Rinse in BSA three times for 5 minutes each.
- Rinse in PBS three times for 5 minutes each.
- Rinse in dH₂O 6 times for 1 minute each.
- Stain with 2 % uranyl acetate (made in dH₂O) for 2 minutes.

TISSUE RUN UP FOR LR WHITE RESIN CHECKLIST

All tissues are to be handled with gloves in the hood

Sample No.: _____

Procedure:

Time:

- | | |
|---|---------|
| 1. 4% paraformaldehyde, 0.1% glutaraldehyde _____ | 4-6 h |
| 2. PBS _____ | 15 min |
| 3. 30% EtOH _____ | 15 min |
| 4. 50% EtOH _____ | 15 min |
| 5. 70% EtOH _____ | 15 min |
| 6. 95% EtOH _____ | 15 min |
| 7. 100% EtOH _____ | 15 min |
| 8. LR White _____ | 1 hour |
| 9. Embedding _____ | |
| 10. Polymerize 60° C _____ | 24 hour |

IMMUNOGOLD LABELING OF LR WHITE EMBEDDED SECTIONS

BSA (PBS containing 1% BSA and 1% Tween 20) used for dilution of block and antibodies

Procedure:

Time:

- | | |
|---|-------|
| 1. Normal goat serum _____ | 1 h |
| 2. primary antibodies _____ | 1 h |
| 3. BSA _____ | 5 min |
| 4. secondary antibodies _____ | 1 h |
| 5. BSA _____ | 5 min |
| 6. PBS _____ | 5 min |
| 7. dH ₂ O _____ | 1 min |
| 8. 2% uranyl acetate in dH ₂ O _____ | 2 min |